




Sentiero... Advanced

The all in one instrument for
OAE, ABR, ASSR and Audiometry



With TE + DPOAE + ABR + ASSR!




Sentiero Advanced

For the first time a hand held, portable lightweight device combines
**Tone and speech audiometry, TEOAE, DPOAE, AEP
(ABR & ASSR)** in a selectable modular format.

The emphasis is on test flexibility, ease of use, reliable diagnostics and minimum test time.

The new **binaural diagnostic ABR** is typically used for NHS conformational diagnostics, neurological diagnostics and objective threshold assessment.



Sentiero... Advanced

The complete solution for your audiology clinic!

Compact & portable handheld, battery powered
Easy to use Colour touch screen, menu guided
Flexibility Configurable test battery: PTA, MAGIC, SPEECH, TEOAE, DPOAE, AEP (ABR & ASSR)
Easy interpretation
Re-editable peak markers on device and PC
Quality parameters
indicating statistical significance (traffic-light)

Detailed test results & test history
up to 1000 tests stored
Data and Patient management Software
MIRA (PC-based)
Easy Interface
USB port for data communication
Interchangeable transducers
Calibration data stored on transducer



Distributed by:

HEALTHWORLD INTERNATIONAL PHILIPPINES CORPORATION

0421 Quirino Avenue La Huerta, Parañaque City, Philippines
Telephone: USA +810 4962007 Philippines +632 8290175 Fax: +632 8266770 email: info@healthwordintl.com