



HEALTHWORLD INTERNATIONAL PHILS, CORP.



EARLAB

■ Hearing aids maintenance station

**EarLab brings out
the professional in you.**



Distributed by: **HEALTHWORLD INTERNATIONAL PHILS, CORP.**

6th Floor Berma Centre, Redemptorist Rd., Baclaran Parañaque City, Philippines 1700
Telephone: USA +810 4962007 PHILIPPINES +632 8792251 Fax: +632 8792111 • email: info@healthworldintl.com



EARLAB

■ Hearing aids maintenance station



Useful therefore indispensable.

All the necessary tools for the revision and cleaning of all type of hearing aids are enclosed in a compact, elegant design:

- the ultrasound bath (fig. 1) and the tunnel to be used as a vacuum bell (fig. 2)
- milliamperometer to control the current usage - 1,2 V max 200 mA - and two syringes, one to aspire and one to blow a fine air stream (fig. 3)
- directional LED light to illuminate your work space with a bright, white light (fig. 3)
- anechoic chamber with speaker (fig. 4)
- stethoscope to acoustically check and compare two hearing aids (fig. 5)

EarLab is the new system for cleaning and checking hearing aids.

EarLab is a compact laboratory on innovative system developed by Labat for the hearing aids dispenser.






EarLab, in tune with your customers.

Give your customers an innovative service by placing their needs on the forefront.

A service that shows your professional know-how and dedication to your customers.

In a few minutes you can change an occasional customer to a loyal repeat customer.



Distributed by:

HEALTHWORLD INTERNATIONAL PHILS, CORP.

The EarLab options.

- Digital microscope with LCD display to check microcircuits (fig. 6)
- Laser Bar Code reader to identify the customer and the hearing aid, to make verbal notes about this appointment and to revive previous contacts (fig. 7)

