




HEALTHWORLD INTERNATIONAL PHILS, CORP.



NASOPHARYNGOSCOPES



Distributed by: **HEALTHWORLD INTERNATIONAL PHILS, CORP.**

6th Floor Berma Centre, Redemptorist Rd., Baclaran Parañaque City, Philippines 1700
Telephone: USA +810 4962007 PHILIPPINES +632 8792251 Fax: +632 8792111 • email: info@healthworldintl.com



NASOPHARYNGOSCOPES



Advertisement addressed to ENT professionals.

See more than others.

The new generation LABAT nasopharyngoscopes convince by their intuitive operability and optimal ergonomics, designed for the tactile sensation and non-tiring handling as required in daily routine operation.

High-End fibre-optic bundles and a high capacity system produce large and bright images at best resolution – necessary for save diagnosis and therapy.

Due to the consequent use and reliable processing of high-class materials, LABAT nasopharyngoscopes are robust and have a long lifetime in spite of the increasing stress of modern cleaning and sterilizing facilities.

- functional design for easy and comfortable handling
- sensitive yet robust control
- optimal flexibility and torsional stiffness of the insertion channel
- high-performance optical system for larger, brighter and sharper images
- suitable reprocessing procedures:
 - fully immersible in cleaning and disinfectant solutions
 - STERIS®
 - ETO
 - machine reprocessing (RDG-E)

Specifications and features can be changed without notice.

	EF-N	EF-NS	EF-N slim	EF-N14	EF-N14 S
diameter of working channel	-	-	-	1.4 mm	1.4 mm
incorporated light-guide cable	-	x	-	-	x
sheath diameter	3.4 mm	3.4 mm	2.8 mm	4.0 mm	4.0 mm
working length	320 mm	320 mm	320 mm	320 mm	320 mm
bending up/down	130°	130°	130°	130°	130°
bending radius	8 mm	8 mm	8 mm	8 mm	8 mm
optical system					
<i>field of view</i>	80°	80°	80°	80°	80°
<i>depth of field</i>	1-50 mm	1-50 mm	1-50 mm	1-50 mm	1-50 mm



Distributed by: **HEALTHWORLD INTERNATIONAL PHILS, CORP.**