

Nasopharyngoscope



Flexible Nasopharyngoscope Model ENT 3.6F

The ENT 3.6F Flexible Nasopharyngoscope is an indispensable device for visualization and diagnosis by otolaryngologists, head & neck surgeons, cosmetic surgeons, speech pathologists, allergists and audiologists.



ENT 3.6FL
with Integral Light Guide



ENT 3.6FP
with Detachable Light Guide

- Small, 3.6 mm diameter shaft for enhanced patient comfort
 - High articulation range for increased maneuverability
- Superior optics and illumination for ample, accurate, color-consistent imaging
 - High-resolution fiberoptic image with a wide field of view
 - Lightweight, sleek, ergonomic design
- Fully immersible, leak testable and video compatible

Specifications:

Flexible Length	30cm
Overall Length	53cm
Shaft Diameter	3.6mm
Articulation Range	135° up / 135° down
Depth of Field	5-50mm
Field of View	70°
Focus Range	+4 to -8 Diopters
Magnification Power	30X



Distributed by: **HEALTHWORLD INTERNATIONAL, INC.**

Telephone: USA +810 4962007 PHILIPPINES +632 8792251 Fax: +632 8792111 • email: info@healthworldintl.com