



# NEW!!!

THE NEWEST, MORE CRISP IMAGE AND SLIMMER DESIGN OF THE PREVIOUS

# BON ENDOCAMERA

Medical Endoscopic Camera  
Very easy to use.....



## Camera Technical Specifications

1. Chip: 1/3" color CCD; 440K pixels
2. Pixels: NTSC: 811(H)×508(V) PAL: 795(H)×595(V)
3. Effective Pixels: NTSC: 768(H)×494(V) PAL: 752(H)×582(V)
4. Sync: Internal Synchronization
5. White Balance: Automatic, Manual, Auto Tracking
6. R and B adjustable
7. Brightness, contrast, sharpness adjustable
8. digital zoom
9. Internal memory which can store up to 4 image pictures
10. Power supply: 100-220V 50HZ
11. Input power: 50VA
12. Fuse: T1A/250V
13. An ordinary video output terminal
14. TV System : NTSC (625 lines; 50Hz)
15. Working condition: +5°C~40°C
16. Humidity: 30%~80%(except the state of dew)
17. Air pressure: 700hPa~1060hPa

## Special Features

SLIMMER  
LIGHT WEIGHT  
BETTER DESIGN  
MORE CRISP IMAGE  
BUILT-IN QUICK CONNECT  
THAT FITS MOST SCOPES  
Electronic Day & Night Control  
Backlight Compensation  
AutoWhite Balance

## Coupler Features / Design

Professional Endoscope Coupler  
Heavy Duty Design  
Fits both Rigid and Flexible Scopes  
High Resolution double coated lens system  
Standard C-Mount and Waterproof ( Fits any C-Mount Cameras)  
Focusing Mechanism  
Spring Loaded Leaf Clamping Mechanism



COUPLER & LENS



CAMERA SYSTEM



SAMPLE IMAGE

Distributed by: **HEALTHWORLD INTERNATIONAL PHILS, CORP.**

6th Floor Berma Centre, Redemptorist Rd., Baclaran, Parañaque City Philippines 1700  
 Telephone: PHILIPPINES +632 8290175 / 3463691 / 09228496753 Fax: +632 8266770  
 email: [info@healthworldintl.com](mailto:info@healthworldintl.com)